

Article	Content
SL-9550-smp	1 mL, 100 rxn á 20µL
SL-9550-5ML	5x 1 mL, 500 rxn á 20µL
SL-9550-10ML	10x 1 mL, 1000 rxn á 20µL
SL-9550-20ML	20x 1 mL, 2000 rxn á 20µL

Storage Conditions**Long-Term Storage**
at -20 °C in the dark**Short-Term Storage**
at 4 °C in the dark

primaQUANT 1STEP MULTIPLEX

2x 1-Step RT-qPCR Probe Master Mix - no ROX

DESCRIPTION

Our **primaQUANT 1STEP MULTIPLEX 2x RT-qPCR Master Mix** is an optimized ready-to-use mixture for probe-based (multiplex) assays such as Taqman®, Beacons and MGBs. It contains all enzymes for both Reverse Transcription and qPCR, as well as dNTPS and MgCl₂ - combined in an optimized buffer system that provides fast kinetics and target amplification even for difficult templates.

The **primaQUANT 1STEP MULTIPLEX 2x Master Mix** contains all components - you just need to add primers, probes and template RNA.

The Master Mix can be used for one-step RT-qPCR, so there is no need for additional reverse transcription.



DID YOU KNOW?

> Some qPCR cyclers require ROX - **primaQUANT 1STEP MULTIPLEX** is also available with low or high concentrations of ROX.

STANDARD PROTOCOL

BEFORE YOU START

- > After thawing, please **invert the Master Mix tube 6-8 times**.
- > **Do not vortex** the Master Mix to prevent damage of the enzyme.

NOTE

- > Cycling conditions highly depend on the primer, probe, amplicon and input material and thus might need adjustments.
- > However, standard cycling conditions can be applied for the majority of qPCR assays out-of-the box.

3-Step qPCR Protocol

Step	Time	Temperature
Reverse Transcription	10-15 minutes	50 °C
Initial Denaturation	1-3 minutes	92 °C - 95 °C
25- 40 cycles		
Denaturation	5 - 10 seconds	92 °C - 95 °C
Annealing	5 seconds	60 °C - depending on primer
Extension	20 - 30 seconds	72 °C

2-Step qPCR Protocol

Step	Time	Temperature
Reverse Transcription	10-15 minutes	50 °C
Initial Denaturation	1-3 minutes	92 °C - 95 °C
25- 40 cycles		
Denaturation	5 seconds	92 °C - 95 °C
Annealing/Extension Combined	20 - 40 seconds	60 °C - depending on primer

ULTRA-FAST PROTOCOL

BEFORE YOU START

- > After thawing, please **invert the Master Mix tube 6-8 times**.
- > **Do not vortex** the Master Mix to prevent damage of the enzyme.

NOTE

- > Ultra-fast Cycling conditions **highly depend on the ramping rate of your qPCR cycler, primer, probe, amplicon and input material** and thus might need adjustments.
- > Ultra-fast cycling conditions can be applied for the majority of qPCR assays out-of-the box, provided that your primer/probe sets do not show unspecific binding.

3-Step qPCR Protocol

Step	Time	Temperature
Reverse Transcription	10-15 minutes	50 °C
Initial Denaturation	60 seconds	92 °C - 95 °C
25- 40 cycles		
Denaturation	1-5 seconds	92 °C - 95 °C
Annealing	1-5 seconds	60 °C - depending on primer
Extension	5 seconds	72 °C

2-Step qPCR Protocol

Step	Time	Temperature
Reverse Transcription	10-15 minutes	50 °C
Initial Denaturation	60 seconds	92 °C - 95 °C
25- 40 cycles		
Denaturation	1 second	92 °C - 95 °C
Annealing/Extension Combined	5-10 seconds	60 °C - depending on primer

RECOMMENDED REACTION MIXTURE PER WELL

Components	20 μ L Reaction	10 μ L Reaction	Final Concentration
2x primaQUANT 1STEP Master Mix	10 μ L	5 μ L	1x
Forward Primer	variable (e.g. 2 μ L)	variable (e.g. 1 μ L)	100 - 600 nM
Reverse Primer	variable (e.g. 2 μ L)	variable (e.g. 1 μ L)	100 - 600 nM
Probe	variable (e.g. 2 μ L)	variable (e.g. 1 μ L)	100 - 600 nM
Template RNA	variable	variable	0.1 - 100 ng/reaction
Sterile Water	adjust to 20 μ L	adjust to 10 μ L	

 NOTE

> For maximum efficiency and specificity annealing temperatures as well as extension time, primer/probe concentration and template concentration need to be optimized.

CALCULATOR TOOL



Please feel free to download our Excel sheet calculator to calculate the necessary volumes:

<http://calculator.steinbrenner-laborsysteme.de>.



qPCR KNOWLEDGE CENTER

<http://www.qpcr-guide.com>

<https://www.steinbrenner-laborsysteme.de>

HOW TO VALIDATE A QPCR SETUP


i BEFORE YOU START

> You can find additional information on how to validate and set up a qPCR in our qPCR Knowledge Center.

Required Controls

- > **RNA Dilution Series**
A RNA dilution series is used to validate the dynamic range, find the optimal RNA input amount and estimate the overall PCR efficiency.
- > **Primer Dilution Series**
High primer amounts can result in unspecific primer binding that limit the fidelity of your qPCR.
- > **No Template Control (NTC)**
A control in which all components except the template are added - this control is used as a negative control and should not show amplification.

Recommended Validation Layout



APPLICATIONS

- Probe-based quantitative PCR
 - TaqMan® Probes
 - Any Dual-Labeled Hydrolysis Probe
 - Molecular Beacons
 - Scorpion Probes
- RNA and miRNA Expression
- Multiplexing (up to 4 colors)
- Transcript Variant Analysis

QUALITY CONTROL

Our **primaQUANT 1STEP MULTIPLEX 2x RT-qPCR Master Mix** undergoes stringent quality controls.

Each lot is tested in a probe-based qPCR with RNA and MS2 Phage RNA input.

Enzyme purity and homogeneity of > 98 % is validated using a Bioanalyzer SDS protein electrophoresis.

All Master Mixes are free of detectable endonuclease- & exonuclease activity:

- Incubation of 1 µg of plasmid DNA with 5 U for 4h at 37 °C and 72 °C
- Incubation of 1 µg of a DNA size standard with 5 U for 4h at 37 °C and 72 °C

FURTHER INFORMATION

For more information, please visit our website

<https://www.steinbrenner-laborsysteme.de>



Steinbrenner
Laborsysteme GmbH

Steinbrenner Laborsysteme GmbH
In der Au 17 · 69257 Wiesenbach · Germany
Fon +49 (0)6223 8612-47
Fax +49 (0)6223 8612-48
mail@steinbrenner-laborsysteme.de
www.steinbrenner-laborsysteme.de



qPCR KNOWLEDGE CENTER

<http://www.qpcr-guide.com>