


Oral sample collection
EASY • PAINLESS • PROVEN

Collect superior samples for your genetic analysis

All expression analysis starts with the collection of RNA samples. ORACollect®-RNA (ORE-100) is an all-in-one system for the collection, stabilization and transportation of RNA from saliva for *in vitro* diagnostic use.

- Quick and reliable collection of RNA
- Easy collection makes it ideal for self or assisted collection
- Stabilizes human, viral and bacterial RNA found in saliva at room temperature
- Compatible with downstream applications (e.g., RT-qPCR, microarray etc.)

Applicable for:

- Human mRNA expression profiling
- Developing “expression signatures” of disease
- Validating detection of infections (e.g., influenza, coronavirus, etc.)

Specifications:

- Catalogue number: ORE-100
- Sample stability: 60 days
- Shelf-life: 12 months


Superior samples • Proven performance

Some DNA Genetek™ products may not be available in all geographic regions.
 *ORACollect is a registered trademark of DNA Genetek Inc. All other brands and names contained herein are the property of their respective owners.
 Patent (www.dnagenetek.com/legalnotices) • © 2020 DNA Genetek Inc., a subsidiary of OraSure Technologies, Inc., all rights reserved. PD-BR-00365 Issue 1/2020-03